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


## TTL/RS-485 datakabel 0.5 M (tillbehör L-OKP050)

# Contents


ADDITIONAL MANUAL FOR CONTROL UNIT MED XW70L	59
Control unit keys	59
Use of temperature control	60
Function display lights	61
Alarms	62
Pr1 user level parameters for temperature and door alarms and their delay times as well as sensor temperature values	63
Copy factory parameters from the Hot Key to the device memory	64
Copy factory parameters from device memory to Hot Key programming stick	65
ADDITIONAL MANUAL FOR CONTROL UNIT MED SAVE XW757K	66
Function display lights	67
Use of temperature control	68
Set the clock time	69
Storage of temperature and alarm data	70
How to export data and alarms to USB	70
Reading temperature and alarm data from USB memory	70
Clearing temperature and alarm information from the devices memory	72
Alarms	73
Pr1 user level parameters for temperature and door alarms and their delay times as well as sensor temperature values	74
Pr2 service and factory level password and parameters	76
Copy factory parameters from the Hot Key to the device memory	77
Copy factory parameters from device memory to Hot Key programming stick	78
Defrost monitor	79
P4 display temperature probe (average temperature)	79
SAVE models - temperature control battery switches	80
DIN 58345 accessory package for MED SAVE control unit	81
Additional NTC probe and RJ45 cabel (client's probe, accessory)	81
Instructions for switching potential-free alarm and testing	82
Retrofitting a potential-free alarm relay (accessory KK9040)	83
TTL/RS-485 data cable 0.5 m (accessory L-OKP050)	84
Notes	85

# ADDITIONAL MANUAL FOR CONTROL UNIT MED XW70L

Single- and double door refrigerators, single door freezer cabinets


Double door combination cabinets

Lower part chiller  
or freezer

Upper part chiller

## Control unit keys

Not in use

Temperature control key

Temperature control key

manual defrost

To display target set point:  
in programming mode it selects a parameter or confirms an operation

LED light

This switch can be used in models with glass doors if having the light on when the door is closed is preferred. Normally, the light switches on and off automatically with the door switch. A signal light displays on the screen when the light has been switched on from the controller switch .

ON/OFF controller

# Use of temperature control

## Switch on the cabinet

- Switch the power on by pressing the key  .

## Adjust temperature

- Press once key  , then the display will show the set temperature.
- Adjust temperature with keys  and  .
- Press again key  or wait for 15 seconds and the display will revert to the normal state.

## Additional defrost (manual defrost)

- Press key  , for more than 3 seconds to start additional defrost.

Additional defrost cannot be terminated by pressing the key. The defrost function stops automatically.

### Note!

**Additional defrost does not turn on if there is no need for defrosting.**

## To lock and unlock the keyboard

- Keep pressed together for more than 3 s the  and  keys. The "POF" message will be displayed and the keyboard will be locked. If the keyboard is locked it is impossible to adjust the temperature or switch off the power of the unit.
- Keep pressed together for more than 3 s the  and  keys till the "Pon" message will be displayed.

## Control of the minimum and maximum temperatures recorded


- Press button  once. "Lo" will appear on the display for a moment followed by the minimum temperature logged.
- Press button  once. "Hi" will appear on the display for a moment followed by the maximum temperature logged.
- The equipment will revert to the normal operation in 5 seconds.

## Reset the the minimum and maximum temperatures recorded

- Press button  "Lo" or  "Hi" and the recorded temperature will appear on the display. Immediately press and hold  button more than 3 s and the "rst" will blink on the display. Recorded temperature will be reset during blinking. The equipment will revert to the normal operation automatically.

Reset both minimum and maximum recorded temperature separately.

## Function display lights

	Light on	Compressor is running
	Light flashing	Time delay to protect compressor after start-up
	Light on	The evaporator fan is on
	Light flashing	Fans delay after defrost in progress.
	Light on	Defrost is on
	Light flashing	Time delay after defrost
	Light on	Internal LED light is on
	Light on	Alarm; see chapter "Alarms"

# Alarms

## Note!

Reset the alarm signal by pressing any button. The alarm code stays on the display.

### **HA High temperature alarm**

The temperature inside the cabinet has exceeded the set alarm temperature. The alarm can be instigated by sudden icing of the evaporator. Initiate additional manual defrost. See instructions on the previous page. Ensure the cabinet has not been loaded with warm products. Check that a door has not been left open. Check that any product is not blocking air circulation within the cabinet. The cabinet is meant to keep ready chilled products cold.

### **LA Low temperature alarm**

The temperature inside the cabinet has fallen below the set value. Ensure that the products put into the cabinet are not too cold.

The temperature alarm will reset automatically after the interior of the cabinet returns to the set temperature range. If this does not happen, move the products to another cabinet and contact your service provider.

### **dA Door open alarm**

Close the door carefully and the alarm will reset automatically.

### **HA2 Condenser alarm (MEDLAB 400 series and freezers only)**

MEDLAB 400 cabinet alarms of dirty condenser / condenser filter. Check that ambient temperature is not more than +32°C and clean condenser / filter. If alarm recurs call to service.

### **P1 Cooling control probe failure \***

### **P2 Defrost control probe failure \***

### **P4 Condensing temperature probe failure (MEDLAB 400 series and freezers only)\***

\*Note! The cabinet operates on its standby system and will hold normal temperature, but reset the alarm signal and immediately contact your service provider.

### **EE Data failure**

To reset an 'EE' Alarm on the display and resume normal function of the equipment press any key for 3 seconds. After this 'rSt' message will be displayed again for 3 seconds and the equipment should return to normal function.

# Pr1 user level parameters for temperature and door alarms and their delay times as well as sensor temperature values

Press simultaneously **SET** + **▽** A text "ALC" will appear on the screen.

To change the parameter ALC, ALU, ALL, ALd or did search with the arrow buttons **△** or **▽** and press **SET** at the desired parameter. Adjust the value with the arrow buttons **△** or **▽** and confirm with selection button **SET**.

The easiest way to turn the alarm on when using factory settings is by leaving the door ajar (door switch activates) for one minute, which results in open door alarm (voltage free alarm relay control)

## Alarms

### ALC alarm mode (rE or Ab)

rE (default factory setting)

- temperature alarm limit changes when adjusting the temperature
  - > for example with the settings Set=4, ALU=4 and ALL=2 the alarm thresholds will be at temperatures of +8°C and +2°C after ALd delay.

### ALU maximum temperature alarm

- setting 4 in both C and F models

### ALL minimum temperature alarm

- setting 2 in C models and 3 in F models

### Ab (option):

- fixed alarm threshold with ALU and ALL remaining unchanged while adjusting temperature
- ALU and ALL are set to preferred alarm threshold

## Delays

- ALd; temperature alarm delay (setting: 20 min)
- did; door alarm delay (setting 1 = 1 minute, 0 = immediate)

## db1-db4 sensor temperature values

Move with the arrow buttons **△** and **▽** to the desired sensor, and press **SET**. The screen will display the temperature of the selected sensor. If the sensor is not in use, a text "noP" appears on the screen. Return to the menu using the button **SET**. The device will return automatically to normal mode within 15 seconds.

## Copy factory parameters from the Hot Key to the device memory

Turn off the machine at the power switch  .

- Display indicates OFF

- Remove the four plastic plugs from the cover, unscrew the screws, and remove the cover.

**Note!**

**Beware of live connectors.**

- Install the programming stick (A) in the Hot Key connector on the controller. The programming stick cannot be inserted the wrong way round.

- Turn on the device with the power switch  .
- The text "doL" flashes on the display, which copies the parameters from the stick to the device memory.

- Finally, the text "End" appears on the display. Allow the text to turn off before removing the programming stick.

- Replace the cover.

**Note!**

If the controller and the programming stick are not compatible, the text "Err" will be displayed.

# Copy factory parameters from device memory to Hot Key programming stick

- The device must be switched on.
- Remove the four plastic plugs from the cover, unscrew the screws, and remove the cover.

**Note!**

Beware of live connectors.

- Install the programming stick (A) in the Hot Key connector on the controller. The programming stick cannot be inserted the wrong way round.

- Briefly press the key to display "UPL". The parameters are copied from the device memory to the programming stick.

- Finally, the text "End" will be displayed. Allow the text to turn off before removing the programming stick.

- Replace the cover.

**Note!**

If the controller and programming stick are not compatible, the text "Err" will be displayed.

# ADDITIONAL MANUAL FOR CONTROL UNIT MED SAVE XW757K

Single- and double door refrigerators, single door freezer cabinets


Double door combination cabinets

Lower part chiller

Upper part chiller

USB connector for copying data stored in the controller

LED light switch (standard in GD models)

Copying the recorded data to a memory stick

Temperature and alarm recording button

Temperature control key  
In programming mode it browses the parameter codes increases the displayed value.

Temperature control key  
in programming mode it browses the parameter codes or decreases the displayed value.

To display target set point:  
in programming mode it selects a parameter or confirm an operation.


Push for 3s to access SET point programming  
Push for 5s to start a manual defrost

ON/OFF controller

## Function display lights

	Light on	Compressor is running
	Light flashing	Time delay to protect compressor after start-up
	Light on	The evaporator fan is on
	Light flashing	Fans delay after defrost in progress
	Light on	Defrost is on
	Light flashing	Drip time in progress
	Light on	An alarm is occurring, see "Larms"
	Light on	Data recording (REC) is on
	Light on	Battery status OK (battery backup)
	Light flashing	Battery is being charged (battery backup)
	Light flashing	Battery is faulty or empty (battery backup)

# Use of temperature control

## 1. Switch on the cabinet

- Switch the power on by pressing the key  3 seconds.

## 2. Adjust temperature

- Checking temperature settings: Press the button  shortly, and the temperature setting displays on the screen.
- Press key  for 3 seconds and "SEt" message will be displayed. Release the button  immediately then the display will show the set temperature.
- Adjust temperature with keys  and .
- Press again key  or wait for 60 seconds and the display will revert to the normal state.

**Note!**

Press key  for more than 5 seconds an additional defrost will be activated (dEF).

## 3. Additional defrost (manual defrost)

- Press the button  for 5 seconds, and the screen will display the text "SET". Keep the button pressed down and the text "dEF" will appear on the screen. The extra thawing will start. The additional thaw function cannot be switched off using the button. Thawing will end automatically.

**Note!**

Additional defrost does not turn on if there is no need for defrosting.

## 4. To lock and unlock the keyboard

- Keep pressed together for more than 3 s the  and  keys. The "POF" message will be displayed and the keyboard will be locked. In lock mode, only minimum and maximum temperatures can be checked. Other functions are locked.
- Keep pressed together for more than 3 s the  and  keys till the "Pon" message will be displayed.

## 5. Control of the minimum and maximum temperatures recorded (P4 additional probe)

- Press button  once. "Hit" will appear on the display for a moment followed by the maximum temperature logged.

- Press button once. "Lot" will appear on the display for a moment followed by the minimum temperature logged.

- wait for 5 seconds and the display will revert to the normal state.

## 6. Reset the the minimum and maximum temperatures recorded

- Press button "Hi" or "Lo" once and the recorded temperature will appear on the display.

Immediately press button more than 3 s and the "rst" will blink on the display.

Recorded temperature will be reset during blinking. The equipment will revert to the normal operation automatically.

Reset both minimum "Lo" and maximum "Hi" recorded temperature separately.

## Set the clock time

Set the clock time by pressing 3 s and message will be displayed

Press and message will be displayed

Confirm by pressing and select HOURS (for example. 9) by buttons och .

Confirm by pressing and go through all the time settings:

Hur = hour

Min = minutes

Udy = day of the week (Sun,Mon,tuE,Ued,thE,Fri,SAt)

dAy = day of the month

Mon = month

yEa = year

Wait for 60 seconds and the display will revert to the normal state.

## Storage of temperature and alarm data

Press the button  for 3 seconds. Use the arrow keys to select  or  parameter "YES" and confirm the function with the button , so the recording starts and the indicator light is displayed . The device stores the information in its memory.


When you want to stop recording, select "no" according to the previous instruction, and the recording stops and the indicator light  the display turns off.

## How to export data and alarms to USB

If you want to copy the data stored on the device, connect the external USB memory device to the USB connector of the device.

**Note!**

2GB external USB-memory (USB 1.1 or 2.0, FAT formatting).

Press the button  for 3 seconds, in which case the stored data is copied to an external USB memory. During copying, texts dAt, Lod, ALr, and Lod are displayed.

When the data is copied, the word "End" is displayed.


If data copying fails, "Err" will be displayed.

## Reading temperature and alarm data from USB memory

Connect the USB memory to your computer and read the USB memory information.

Two files, ALARM and DATA, are stored in USB memory.

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|**IDX, LABEL, DESCRIPTION, START, STOP**

```

1,PWF,Power Failure,12-30-21 11:54,12-30-21 11:54
2,dA,Door Open Alarm,12-30-21 11:22,12-30-21 11:23
3,HA4,High Temperature Alarm Probe TL2,12-30-21 07:58,12-30-21 08:31
4,HA1,TR3 High Alarm,12-30-21 07:58,12-30-21 08:36

```

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IDX	Measurement cycle
DATE	Date, time
TR3	Cooling control probe
EVP	Defrost control probe, evaporator probe
TL1	Condensing temperature control probe
TL2	Display temperature monitoring probe
ST	Compressor status
-	Compressor status: compressor is not running
C	Compressor status: compressor is running
D	Compressor status: defrosting is in progress
OFF	Compressor status: compressor is switched off

# Clearing temperature and alarm information from the devices memory

Press simultaneously  +  for 3 seconds, whereby text "Pr1" will be displayed.

Confirm the selection with  and select with arrows  and  "rSd" or stored temperature data.

Press the button  and use the arrow keys to select  and  "YES". Confirm the selection with the button , whereby the temperature data is erased from the memory.

When the temperature data is erased from the memory, "rSA" i.e., the stored alarms is displayed.

Press the button  and use the arrow button to select  "YES". Confirm the selection with , whereby the alarm information is erased from the memory.

Wait 60 seconds after which the unit will return to normal mode.

# Alarms

## Note!

Reset the alarm signal by pressing any button. The alarm code stays on the display.

### **HA1 High temperature alarm (internal temperature probe 1, TR3)**

The temperature inside the cabinet has exceeded the set alarm temperature. The alarm can be instigated by sudden icing of the evaporator. Initiate additional manual defrost. See instructions on the previous page. Ensure the cabinet has not been loaded with warm products. Check that a door has not been left open. Check that any product is not blocking air circulation within the cabinet. Temperature alarm will be cancelled automatically when the inside temperature has returned to set limits. The cabinet is meant to keep ready chilled products cold.

### **LA1 Low temperature alarm (internal temperature probe 1, TR3)**

The temperature inside the cabinet has fallen below the set value. Ensure that the products put into the cabinet are not too cold.

The temperature alarm will reset automatically after the interior of the cabinet returns to the set temperature range. If this does not happen, move the products to another cabinet and contact your service provider.

### **PFI Cooling control probe failure (internal temperature probe 1, TR3)**

### **PF2 Defrost control probe failure (probe 2)**

### **HA3 Condenser high temperature alarm (condenser probe 3, TL1, MEDLAB 400 series and freezers only\*)**

Cabinet alarms of dirty condenser / condenser filter. Check that ambient temperature is not more than +32°C and clean condenser / filter. If alarm recurs call to service.

### **LA3 Condenser low temperature alarm (condenser probe 3, TL1\*)**

### **PF3 Condensing temperature probe failure (condenser probe 3, TL1\*)**

### **HA4 Display temperature high temperature alarm (additional probe 1, TL2)**

### **LA4 Display temperature low temperature alarm (additional probe 1, TL2)**

### **PF4 Display temperature probe failure (additional probe 1, TL2)**

### **dA Door open alarm**

Close the door carefully and the alarm will reset automatically.

# Pr1 user level parameters for temperature and door alarms and their delay times as well as sensor temperature values

## 1. Pr1 user level parameters for temperature and door alarms

Press simultaneously  +  buttons for about 3 seconds whereby text "Pr1" will be displayed.

Press the button  and use the arrow keys  and  to search parameter itP, AIC, AIU, AIL, Ald, A4U, A4L, A4d or did to be changed.

Press the button  at the parameter to be changed.

Change the value with the arrow button  and  . Confirm your selection with  button.

By default, the alarm can be accessed by holding the door ajar (door switch releases) for one minute, after which the door-open alarm releases (checking for potential-free alarm relay).

**itP = recording interval**

**AIC = alarm mode**

**AIU = upper limit for alarm P1**

**AIL = lower limit for alarm P1**

**AIL = alarm delay P1**

**A4U= upper limit for alarm P4**

**A4L = lower limit for alarm P4**

**A4d = alarm delay P4**

**did = door alarm delay**

## 2. P1 cooling control probe Alarms

### AIC alarm mode (rEL or AbS)

**rEL** (selected as factory default):

- the temperature alarm limit changes when the temperature changes.

For example, with the setting Set = 4 and AIU = 4 and AIL = 2, alarms come at + 8C and + 2 ° C after Ald delay

### AIU maximum temperature alarm

- setting value 4 for C and F models

### AIL minimum temperature alarm

- setting value 2 for C models and 3 for F models

### AbS (selectable):

- fixed alarm limit, whereby AIU and AIL do not change when the temperature is changed
- AIU and AIL are set to preferred alarm threshold

## 3. P4 Display temperature probe alarms (average temperature)

P4 sensor has always set (AbS) alarm thresholds, which are set with parameters A4L and A4U:

Factory value:	MED chiller		MED freezer	
	A4L	A4U	A4L	A4U
	1	19	-30	-13

## 4. Delays

- Ald temperature alarm delay (setting value 20 min) P1
- A4d temperature alarm delay (setting value 10 min) P4
- did Door alarm delay alarm (Setpoint 1 = 1 minute, 0 = immediately)

## 5. db1-db4 temperature values

**db1 = temperature control probe P1**

**db2 = defrost control probe (vapor probe) P2**

**db3 = condensing temperature probe /condenser filter P3**

**db4 = display temperature monitoring (average temperature, so-called P4 monitoring sensor)**

Press the button  for 3 seconds and use the arrow keys  and  to move to the desired sensor and press  , whereby the temperature value of the selected sensor is displayed.

Press the button  , to access the next parameter. If the sensor is not in use, text "nPr" is displayed.

Keep browsing the menu by  button. The unit will automatically return to normal mode within 60 seconds.

## Pr2 service and factory level password and parameters

### Pr2 service level password 012

Press simultaneously  and  buttons for about 3 seconds whereby text "Pr1" will be displayed.

Press the button , whereby text "Pr2" will be displayed.

Press  and enter the password. 0-- flashes on the display.

Press  and select number 1 by pressing the key , whereby the display will show.

Press  and press the key  to select number 2, whereby 012 will be displayed.

Press the button , and you are at level Pr2.

### Service and factory level Pr2 parameters

Search the parameter to be changed by arrow keys  and .

Press  for the parameter to be changed.

Change the value with the arrow keys  and  . Confirm the selection by  button

The unit will automatically return to normal mode within 60 seconds.

### Device parameters:

Festivo Service & Spare Parts Sales  
Tel. +358 3 553 8681  
service@festivo.fi

or an authorized service center

## Copy factory parameters from the Hot Key to the device memory

- The device must be switched on.
- Remove the four plastic plugs from the cover, unscrew the screws, and remove the cover.

**Note!**

**Beware of live connectors.**

- Install the programming stick (A) in the Hot Key connector on the controller. The programming stick cannot be inserted the wrong way round.

- Press the key for three seconds to display "HA".

- Use the arrow keys and to find the parameter "doL".

- Press the button until "doL" starts flashing. The parameters are copied from the stick to the device memory. Finally, the text "End" appears on the display. Allow the text to turn off before removing the programming stick.

- Replace the cover.

**Note!**

If the controller and programming stick are not compatible, "Err" will be displayed.

# Copy factory parameters from device memory to Hot Key programming stick

- The device must be switched on.
- Remove the four plastic plugs from the cover, unscrew the screws, and remove the cover.
- Note!** Beware of live connectors.
- Install the programming stick (A) in the Hot Key connector on the controller. The programming stick cannot be inserted the wrong way round.

- Press the key for three seconds to display "HA".

- Use the arrow keys and to find the parameter "UPL".

- Press the button until the UPL starts to flash. The parameters are copied from the device memory to the programming stick. Finally, the text "End" appears on the display. Allow the text to turn off before removing the programming stick.

- When "End" appears on the display, you can remove the Hot Key stick.
- Replace the cover.

**Note!**

If the controller and programming stick are not compatible, the text "Err" will be displayed.

# Defrost monitor

## Defrost monitor operation

Festivo MEDLAB products have automatic defrosting. Defrosting is controlled by the cabinet's own monitoring sensor P2 and defrosting is stopped as soon as the evaporator has melted. The defrost monitor function is controlled by the bottle sensor P4.

The defrost monitor is a function that protects the storage of highly sensitive products (e.g., reagents) during the defrost cycle. The defrost monitor monitors the internal temperature during the defrost cycle and stops defrosting if the internal temperature rises above the set limit, for example, due to the door opening.

If the defrost monitor function interrupts defrosting, a new defrost cycle will start automatically after 1 hour when the temperature has stabilized.

Defrost monitor factory settings:

Medicine refrigerators (+ 7°C)

Freezers (-10°C)

## Defrost monitor and parameter Srt setting (Pr2 service level)

The defrost monitor is set at pr2 level with parameter Srt.

### Note

**The value of the Srt parameter does not change even if the cabinet setpoint (SET) is changed.**

### Note!

**The value of the parameter Srt must not be too close to the set (SET) temperature.**

Medicine refrigerators: set temperature (SET) + 2 or more = value of parameter Srt

Freezers: set temperature (SET) + 4 or more = value of parameter Srt

For example, temperature setpoint SET = -22

-22 + 4 = -18 (Srt)


## P4 display temperature monitoring probe (average temperature)

The P4 control sensor is installed in the bottle on the rear wall of the cabinet. Loosen the bottle cap to remove it from the rack.

## SAVE models - temperature control battery switches


MEDLAB SAVE 400 and MED SAVE standard models with a press switch

A1 = 1-door devices and  
2-doors lower cabinet

A2 = 2-doors upper cabinet

MEDLAB SAVE 400 DIN accessory and BB standard model with lockable key switch

B1 = 1-door devices and  
2-doors lower cabinet

B2 = 2-doors upper cabinet

### Standard equipment:

Switch on the battery from the black push switch (A) on the back panel of the appliance. The controller will alert you if the battery is not connected.

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
Switch on the battery from the lockable key switch (B). Turn the key from the back to the right. Do not leave the key in lock, but keep the key in the agreed location. Battery backup can be turned off during transport or transfer of the unit. Turn the key 90° to the left, so that battery backup is switched off.

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There is a potential-free connection point (A) at the back of the cabinet for alarms. DIN 58345 equipment package includes a battery backup for MED SAVE controller power failure and temperature registering during 12 hours, lockable battery back-up user switch (B), and antifreeze function.

BB cables have battery back-up by default.

## Additional NTC probe and RJ45 cable (client's probe, accessory)

- Remove the protective cap (E)
- Twist to open the bushing of the push through plug (F)
- Push sensor (G) through the push through hole into the box
- Plug in the cables into the switch strip (C)
- Twist to close the push through plug (F)

### Note!

Customer sensor must be calibrated according to the box temperatures. Sensor calibration depends on its location.

## RJ45 connection cable

### Alarm switch factory setting:


4/blue

5/blue-white

Also a separate sensor for continuous temperature control can be connected to RJ45 cable, which can then be monitored through the property control centre (VAK).

We recommend a customer sensor installation at the factory.

Connect RJ45 cable (D) to connector A (connector A, see picture on previous page).

Cable colour codes:

		Pin		Pin	RJ45	
	orange/white	1		1	orange/white	
	orange	2		2	orange	
	green/white	3		3	green/white	
	blue	4		4	blue	
	blue/white	5		5	blue/white	
	green	6		6	green	
	brown/white	7		7	brown/white	
	brown	8		8	brown	

## Instructions for switching potential-free alarm and testing

For FESTIVO MED models (XW70L), potential-free switching readiness is an accessory. For MEDIN (XW777K) and MED SAVE (XW757K) models, the switching readiness is standard.

Behind the product is a potential-free switch bar (connector A, picture on the previous page):

NC = Normal closed

NO = Normal open

C = Common

XW70L =AOP

XW777K =AOP, Aro

XW757K =AOP, Aro

The effect of a power failure on a potential-free alarm can be set in Pr2 level parameters with the parameters AOP and Aro. The AOP parameter selects the operation of the alarm relay (closing/opening). The Aro parameter is used to turn the battery/alarm relay function of a power failure on/off.

Parameter AOP

oP=Open cL= Close (XW70L)


oP=Open cL= Close( XW757K , XW777K ) , Battery Back-up (XW757K , XW777K ) , Aro = YES

For battery back-upped devices, parameter Aro=YES so that the alarm also works during a power failure.

**Use the multimeter to test the top data from connector A.**

## Retrofitting a potential-free alarm relay (accessory KK9040)


1. Remove the cabinet cover plate. Four fastening screws are under white protective plugs.
2. Install the relay/relays (1) and the terminal strip (2) in the position shown below with velcro and connect the relay to the terminal strips.

One control unit

Two control units

Terminal strip NO NC C

## ONE CONTROL UNIT

K1

## TWO CONTROL UNITS

K1

K2

11-12 NC

11-14 NO

21-22 NC

21-24 NO

TTL/RS-485 data cable 0.5 m (accessory L-OKP050)

## Notes

## Notes

## Notes

Festivo appliances are manufactured in Hollola, Finland, where our research & development, factory, administration and service are also situated.

Festivo is a genuine Finnish product. The Festivo label on the front panel of your cooling appliance is guarantee of high Finnish quality, evolving from long experience, strong materials, and construction, proven over decades.

**Suomen Kotikylmiö Oy**  
Ravitie 3, FIN-15860 Hollola  
Tel. +358 3 553 8600  
[festivo@festivo.fi](mailto:festivo@festivo.fi)

**Festivo Service**  
Keskikankaantie 17, FIN-15860 Hollola  
Tel. +358 3 553 8681  
[service@festivo.fi](mailto:service@festivo.fi)